



14425



14426 (Manual Pump)



14427 (Electric Pump)

COMPARE TO:  
Drive Tender Cloud®  
14426 AQTC-8035\*

**250 LB**  
WEIGHT CAPACITY

### Premium Guard - Deluxe Static Air Mattress Overlay

14425 35" (W) x 72" (L) x 2.5" (H),  
Static Air Mattress, 1/cs  
14426 Manual Pump for Air Mattress, 1/cs  
14427 Electric Pump for Air Mattress, 1/cs

**HCPCS code: E0181**

#### FEATURES:

- Over 240 air cells and chambers "self-adjust" for pressure equalization.
- Three layer system adjusts to evenly distribute air volume.
- Heavy duty gauge vinyl provides maximum air loss protection.
- Easy to clean with mild detergent.
- Four heavy-duty, elastic bed straps prevent mattress from moving.
- Emergency air release valve deflates mattress quickly and safely.
- Compatible with most air filling devices.
- 90 Day Limited Warranty on seams and valve.

#### SPECIFICATIONS

Item #	Width	Depth	Height	Product Weight	Carton Shipping Weight
14425	35"	72"	2.5"	8 lbs.	11 lbs.

• Not endorsed or authorized by Graham Field®.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.