

250 LB

WEIGHT CAPACITY



14428

COMPARE TO:

drive EHOB®

14428 1006EC*



14427 (Electric Pump)



14426 (Manual Pump)




14400

250 LB

WEIGHT CAPACITY

Fill Hose



NEW PACKAGING!



14401

Premium Guard - Static Air Mattress Overlay

- 14428 35" (W) x 72" (L) x 2.5" (H), Static Air Mattress, 1/cs
- 14426 Manual Pump for Premium Guard Air Mattress, 1/cs
- 14427 Electric Pump for Premium Guard Air Mattress, 1/cs

HCPCS code: E0181

FEATURES:

- Unique construction allows for efficient air flow throughout mattress for use in the treatment and prevention of pressure sores.
- Easy to clean with mild detergent.
- Emergency air release valve deflates mattress quickly and safely.
- Air flow holes reduce heat and fluid build up.
- Adjustable bed straps prevent mattress from moving.
- Two adapters for easy filling and air removal.
- 90 Day Limited Warranty on seams and valve.

Premium Guard - Water Mattress

- 14400 32" (W) x 68" (L) x 4" (H), Water Mattress, 1/cs
- 14401 Fill hose for Water Mattress (Sold separately)

FEATURES:

- Three water baffles enhance movement and support to reduce pressure and conform to varying body weights.
- Filling lines and "Foot Placement" lines provide for proper patient positioning.
- Double sealed seams provide maximum strength and durability.
- Vinyl construction is durable and easy to clean with mild detergent.
- Fill lines are conveniently located.
- Size: 68" x 32"
- Limited 90 Day Warranty on seams and valve.

SPECIFICATIONS

Item #	Width	Depth	Height	Product Weight	Carton Shipping Weight
14428	35"	72"	2.5"	8 lbs.	11 lbs.
14400	32"	68"	4"	4 lbs.	6 lbs.

* Not endorsed or authorized by Graham Field®.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

pressure prevention

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM