

## Premium Guard Gel/Foam Overlay

14893	Gel Overlay, 76" (L) x 34" (W) x 3.5" (H), 275 Wt. Cap., 1/cs
HCPCS code: E0185	
14901	Bariatric Gel Overlay, 76" (L) x 42" (W) x 3.5" (H), 500 Wt. Cap., 1/cs
HCPCS code: E0185	

### FEATURES:


- Ideal combination of 3 gel bladders, high density foam shell and support base.
- Optimal prevention, treatment and management of pressure ulcers.
- One piece construction requires no assembly.
- 4 Elastic corner straps secures overlay to mattress.
- Urethane coated nylon cover is, water resistant and vapor permeable.
- Vinyl base is stain resistant and waterproof.
- Bladders contain an aqueous non-toxic gel.
- Designated head and foot sections.
- 1 Year Limited Warranty on cover and mattress.

**SURFACE COMFORT  
TOP LAYER**

**GEL LAYER**

**GEL OVERLAY**

**PREMIUM GUARD OVERLAY  
WITH SOFT NYLON COVER!**



14904

## Bed Side Rail Safety Pads

14904 70" (W) x 13" (D) x 1" (H), 1pr/cs

### FEATURES:

- Prevents semi-conscious, comatose and overdosed patients from injuring their limbs against the bed side rails.
- Manufactured with foam core.
- Covered in vinyl with velcro straps to securely fasten to bed side rails.
- Folds in three sections for convenience of use and storage.

### SPECIFICATIONS

Item #	Width	Depth	Height	Product Weight	Carton Shipping Weight	Weight Limit
14893	34"	76"	3.5"	47.8 lbs.	50.8 lbs.	275 lbs.
14901	42"	76"	3.5"	50.5 lbs.	53.5 lbs.	500 lbs.
14904	70"	13"	1"	2.3 lbs.	3.3 lbs.	NA

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

pressure prevention

To order call toll free: 877.224.0946

**drive**  
WWW.DRIVEMEDICAL.COM