


**COMPARE TO
JAY® BASIC**



14889

General Use Back Cushion with Lumbar Support

14889 18" (W) x 17" (D), 300 Weight Cap., 1/cs

14906 16" (W) x 17" (D), 300 Weight Cap., 1/cs

HCPCS code: E2611

FEATURES:

- Improves seating posture while providing back support and comfort.
- Conforms to provide support to the lumbar region.
- Sealed liner protects foam core from moisture.
- The stabilization board secured on back of foam core prevents hammocking in sling-back wheelchairs.
- Has a non-slip vinyl back to prevent the cushion from sliding.
- 18 Month Limited Warranty on cover and cushion.

**COMPARE TO
INVACARE®
CAREGUARD 101**



15019

Therapeutic 5 Zone Support Mattress

15019 35" (W) x 80" (D) x 5.5" (H), (Has replaced 15019CA), 300 Weight Cap., 1/cs

HCPCS code: E0184

FEATURES:

- Deluxe foam mattress with convoluted top layer reduces pressure.
- High resilient firm foam base to provide support and comfort.
- Deluxe nylon cover is durable, zippered, waterproof anti-microbial and meets California Flame Retardant Standards.
- Deluxe nylon cover reduces friction and shear.
- 2 Year Limited Warranty on foam.
12 Month Limited Warranty on cover.

SPECIFICATIONS

Item #	Width	Depth	Height	Product Weight	Carton Shipping Weight	Weight Limit
14889	18"	17"	2.5"	2.1 lbs.	3.0 lbs.	300 lbs.
14906	16"	17"	2.5"	2 lbs.	3.0 lbs.	300 lbs.
15019	35"	80"	5.5"	16 lbs.	19 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.