



**NEW
IMPROVED
ARM!**

Steel Drop-Arm Commode with Wheels and Padded Arms

- 11101W-2 Assembled, with Wheels, 2/cs
- 11101WL Caster without Lock for 11101W-2, 2/bx
- 11101WLL Caster with Lock for 11101W-2, 2/bx

FEATURES:

- Easy-to-release arm mechanism allows for safe lateral patient transfers to and from commode.
- Arm provides additional support while transferring.
- Manufactured with 1" steel.
- White powder coated frame is attractive and easy to maintain.
- Padded arm rests provide extra comfort and support.
- Comes with 7.5 qt. commode bucket with carry handle and splash shield.
- Rust resistant, swivel casters (2 with locks)
- **Limited Lifetime Warranty.**



WHEN ARM IS IN FOLDED POSITION, THE ARM FRAME PROVIDES A "HANDLE" SUPPORT FOR TRANSFERS.



NEW REINFORCED TOOL FREE BACK REST DESIGN.



**NEW
IMPROVED
CASTERS!**



NON-SKID, RUST RESISTANT, SWIVEL CASTERS. TWO WITH LOCKS!



ENVIRONMENT FRIENDLY PRODUCT
Recycled materials were used to make this product.

Aluminum Commode with Wheels

- 11112-1* Assembled, 1/cs
- 11112-2 Assembled, 2/cs

FEATURES:

- 3" Non-skid, rust resistant, swivel casters (2 with locks).
- Anodized aluminum frame.
- Back easily removes, tool free.
- Angled back for support and comfort.
- Full length plastic armrests.
- Comes with a 7.5 qt. commode bucket, splash shield and cover. Splash shield acts as pail holder.
- **Limited Lifetime Warranty.**

SPECIFICATIONS

Item #	Width (between arms)	Outside Width	Width Seat	Depth Seat	Height (highest)	Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
11101W-2	18"	24"	14"	16.5"	24"	18"	22.5 lbs.	48 lbs.	300 lbs.
11112-1/2	18.5"	24"	14"	16.5"	18.75"	18.75"	9 lbs.	11/21 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.