



ENVIRONMENT FRIENDLY PRODUCT

Recycled materials were used to make this product.



11125KD-1



TOOL FREE ASSEMBLY



11125PSKD-1



TOOL FREE ASSEMBLY



11125PSKD-1

SHOWN WITH DROPPED DOWN ARM

K.D. Deluxe Steel Drop-Arm Commode (Tool Free)

- 11125KD-1 New Knock-Down Frame, (Has replaced 11125-1), 1/cs
- 11125KD-2 New Knock-down Frame, (Has replaced 11125-2), 2/cs
- 10386G-8 Replacement Tip, 1 1/8" Grey, (Has replaced 11123G-8), 8 bx/cs

FEATURES:

- Easy to assemble frame.
- Easy-to-release drop arm mechanism allows for safe lateral patient transfers to and from commode.
- Ideal for those with limited dexterity.
- Legs are height adjustable.
- Easy to clean grey powder coated steel finish.
- Padded armrests and back provide extra comfort and support.
- Comes with 12 qt. commode bucket with carry handle and splash shield.
- **Limited Lifetime Warranty.**


K.D. Deluxe Steel Drop-Arm Commode with Padded Seat (Tool Free)

- 11125PSKD-1 New Knock-Down Frame, with Padded Seat with Open Front End, 1/cs
- 10386G-8 Replacement Tip, 1 1/8" Grey, (Has replaced 11123G-8), 8 bx/cs
- 11125-PS Replacement Seat, 1/cs

FEATURES:

- Easy to assemble frame.
- Padded open front vinyl toilet seat.
- Easy-to-release drop arm mechanism allows for safe lateral patient transfers to and from commode.
- Ideal for those with limited dexterity.
- Legs are height adjustable.
- Easy to clean grey powder coated steel finish.
- Padded armrests and back provide extra comfort and support.
- Comes with 12 qt. commode bucket with carry handle and splash shield.
- **Limited Lifetime Warranty.**

TOOL FREE ASSEMBLY IN 4 EASY STEPS!



SPECIFICATIONS

Item #	Width (between arms)	Outside Width	Width Seat	Depth Seat	Height (highest)	Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
11125KD-1/2	17.5"	24"	13.5"	16.5"	23"	19"	18.1 lbs.	20.1/39.2 lbs.	300 lbs.
11125PSKD-1	17.5"	24"	13.5"	16.5"	23"	19"	20.2 lbs.	22.2 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

commodes

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM