

REFER TO PAGE 182-190
FOR LIFT ACCESSORIES
AND SLINGS.

RECOMMENDED
FOR USE WITH 2 &
4 POINT SLINGS!

Battery Powered Patient Lift with Six Point Cradle

13240 1/cs

BARIATRIC
450 LB
WEIGHT CAPACITY

FEATURES:

- Fast, easy charge from AC outlet.
- **To charge, simply plug in unit. Battery does not need to be removed to be charged. Lift will not operate while plugged in.**
- Heavy gauge construction.
- Jumbo actuator pump with emergency button can switch to manual mode lowers patient down safely.
- Audible warning when battery power falls below 50%.
- The "CHARGE" light will automatically shut off when battery is fully charged.
- Foldable mast and base for easy storage and transportation.
- **Ideal for Homecare use.**
- For use with 4 point and 2 point slings.
- Wide base provides additional stability.
- Lifting capacity: 450 lbs.
- Limited Lifetime Warranty.



NEW & IMPROVED
CRADLE

13240
with Sling
(Sold Separately)

CAN FIT UNDER
A LOW BED!

Battery Powered Patient Lift with Rechargeable, Removable Battery

13242 1/cs

BARIATRIC
450 LB
WEIGHT CAPACITY

FEATURES:

- Same features as 13240.
- Includes removable, rechargeable battery, (item #13258) external charger (item #13255) and wall mounted bracket (item #13260).
- **Ideal for institutional use.**

To avoid lift downtime - order an additional re-chargeable battery, item # 13258. One for the lift and one for the charging station.

REPLACEMENT PARTS:



**External Charger
with Power Cord**

13255 1/ea



**Wall Mounted
Bracket**

13260 1/ea



**Rechargeable
Battery**

13258 1/ea



NEW & IMPROVED
CRADLE

13242
with Sling (Sold Separately)

CAN FIT UNDER A LOW BED!

SPECIFICATIONS

Item #	Boom Maximum Height	Boom Minimum Height	Base Width (open)	Base Width (closed)	Base Height	Base Clearance	Casters	Product Weight	Carton Shipping Weight	Weight Limit
13240	76"	27.5"	41"	15"	4.5"	2"	3" Front/4" Rear	106 lbs.	116 lbs.	450 lbs.
13242	76"	27.5"	41"	15"	4.5"	2"	3" Front/4" Rear	106 lbs.	116 lbs.	450 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

patient room

drive
WWW.DRIVEMEDICAL.COM

To order call toll free: 877.224.0946