



**NEW & IMPROVED CRADLE**

13244  
with Sling  
(Sold Separately)

**CAN FIT UNDER A LOW BED!**

**BARIATRIC**  
**600 LB**  
WEIGHT CAPACITY

**REFER TO PAGE 182-190 FOR LIFT ACCESSORIES AND SLINGS.**

**Bariatric Battery Powered Lift with Six Point Cradle**

13244 1/cs

**FEATURES:**

- 600 lb. weight capacity.
- To charge, simply plug in unit. Battery does not need to be removed to be charged. Lift will not operate while plugged in.
- Jumbo actuator pump with emergency button can switch to manual mode lowers patient down safely.
- Base legs adjust easily and lock securely into open position with padded shifter handle.
- **Ideal for homecare use.**
- Extra wide coated cradle with 360° rotation with 6 point hook up.
- High lift range allows transfer to high surfaces.
- Dual, 3" heavy duty front casters.
- Warranty: Lifetime on frame, 1 year on components and electrical parts.

**NEW & IMPROVED CRADLE**

**BARIATRIC**  
**600 LB**  
WEIGHT CAPACITY



**CAN FIT UNDER A LOW BED!**

13245  
with Sling  
(Sold separately)

**Bariatric Battery Powered Patient Lift with Rechargeable, Removable Battery**

13245 1/cs

**FEATURES:**

- Same features as 13244
- Includes removable, rechargeable battery, (item #13258) external charger (item #13255) and wall mounted bracket (item #13260).
- **Ideal for institutional use.**

*To avoid lift downtime - order an additional re-chargeable battery, item # 13258. One for the lift and one for the charging station.*

**REPLACEMENT PARTS:**



**External Charger with Power Cord**

13255 1/ea



**Wall Mounted Bracket**

13260 1/ea



**Rechargeable Battery**

13258 1/ea

*\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.*

**SPECIFICATIONS**

Item #	Boom Maximum Height	Boom Minimum Height	Base Width (open)	Base Width (closed)	Base Height	Base Clearance	Casters	Product Weight	Carton Shipping Weight	Weight Limit
13244	74.5"	32"	45.75"	25.75"	4.25"	2"	3" Front/4" Rear	127 lbs.	139 lbs.	600 lbs.
13245	74.5"	32"	45.75"	25.75"	4.25"	2"	3" Front/4" Rear	127 lbs.	139 lbs.	600 lbs.