


# SMALL BASE QUAD CANES



## Quad Canes, Small Base

10301-4	Vinyl Contoured Grip, Chrome Cane, 4/cs
10301F-4	Foam Rubber Grip, Chrome Cane, 4/cs
10310-4	Vinyl Contoured Grip, Black Cane, (Has replaced 10310F-4), 4/cs
10311-4	Vinyl Contoured Grip, Bronze Cane, (Has replaced 10311F-4), 4/cs
10311F-4*	Foam Rubber Grip, Bronze Cane, 4/cs

## Quad Canes, Small Base, Ortho Grip

10301G-4*	Ortho K Grip, Chrome Cane, 4/cs
-----------	---------------------------------

### FEATURES:

- 4 point base, combined with offset handle provides additional stability and support.
- Easy-to-use, one-button height adjustment with locking bolt, allows quick, safe fitting.
- Height adjustment 30" - 39".
- No assembly required.
- 300 lb. Weight Capacity.
- Limited Lifetime Warranty.

