

BARIATRIC
400 LB
WEIGHT CAPACITY

drive
1 THIS ITEM IS AVAILABLE IN A 1-PACK!

EASY HEIGHT ADJUSTMENT WITH SELF THREADING KNOB!




10215RD-1
SEE PAGE 52 & 53 FOR REPLACEMENT SEATS.

Go-Lite Bariatric Steel Rollator, Padded Seat, 8" Casters with Loop Locks

10215BL-1 Blue, (Has replaced 685 series and 466), 1/cs
10215RD-1 Red, (Has replaced 685 series), 1/cs

FEATURES:

- Large 8" casters are ideal for indoor and outdoor use.
- Comes with large 21" (W) x 9" (D) x 6" (H) basket that can be mounted under seat (item # 10215BSK).
- Tool Free removable padded back rest for comfort.
- Strong Steel reinforced frame supports 400 lbs.
- Soft padded oversized seat (18" (W) x 13" (D)).
- Special loop lock made of internal aluminum casting operates easily and ensures safety.
- Handles are adjustable in height.
- Limited Lifetime Warranty. Lifetime Warranty on brake cable.



PARTS DIAGRAMS ON PAGE 67

Adjustable Width Winnie Walker, Bariatric EXP Steel Rollator, 8" Casters with Loop Locks

466* Light Blue & Chrome, 1/cs

FEATURES:

- Large 8" casters are ideal for indoor and outdoor use.
- 4" adjustable width allows entry into narrow doorways.
- Strong Steel reinforced frame supports 400 lbs.
- Special loop lock made of internal aluminum casting operates easily and ensures safety.
- Soft padded oversized seat.
- Handles are adjustable in height.
- Comes with oversized 14" x 12" x 7" wire basket (item # 828).
- Limited Lifetime Warranty. Lifetime Warranty on brake cable.

BARIATRIC
400 LB
WEIGHT CAPACITY

EASY TO REACH KNOB ADJUSTS ROLLATOR WIDTH 4" TO ALLOW ENTRY INTO NARROW DOORWAYS.



466*
SEE PAGE 52 & 53 FOR REPLACEMENT SEATS.

ADJUST FRAME WIDTH UP TO 4" WHEN USING LEVER



**SEE ROLLATOR ACCESSORIES ON PAGES 50-57.
ROLLATOR PARTS DIAGRAMS ON PAGE 58-67.
ROLLATOR COLOR PAINT SWATCHES ON PAGE 68.**

SPECIFICATIONS

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

Item #	Handle Height (highest)	Handle Height (lowest)	Unit Width	Unit Length	Seat (depth)	Seat (width)	Seat (height)	Product Weight	Carton Shipping Weight	Weight Limit
GO-LITE BARIATRIC	39"	34"	30"	29"	13"	18"	22"	26 lbs.	30.5 lbs.	400 lbs.
EXP/466	38"	34"	25" - 29"	30"	10"	20"	25"	28 lbs.	32.5 lbs.	400 lbs.

rollators

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM