



WEIGHTS ONLY 14 LBS.!



PARTS DIAGRAMS ON PAGE 59 & 60

Aluminum Rollator with Fold Up and Removable Back Support, Removable 8" Wheels, Padded Seat and Loop Locks

728RD-RTL Red, Retail, 1/cs

FEATURES:

- Comes with new seamless padded seat (14" x 14") with zippered pouch under seat.
- 8" Black non-marring removable casters, ideal for indoor and outdoor use.
- Removable hinged, padded backrest can be folded up or down.
- Comes with standard carry pouch under the seat for added privacy and security of personal items.
- Easy to use deluxe loop locks.
- Ergonomic handles are easy to grip and are adjustable in height.
- Limited Lifetime Warranty.



COMES WITH NEW SEAMLESS PADDED SEAT WITH CARRY POUCH!

OPTIONAL BASKET AVAILABLE MODEL # 9501S72606

SEE PAGE 52 & 53 FOR REPLACEMENT SEATS.

EASY HEIGHT ADJUSTMENT WITH SELF THREADING KNOB!



CERATED BRAKES



ERGONOMIC HANDLE FOR EASY GRIP!



FOLD-UP & REMOVABLE BACKREST!



728RD-RTL



PADDED SEAT OPENS TO A ROOMY CONVENIENT STORAGE POUCH! ALSO COMES WITH INSIDE ZIPPER POUCH ON UNDERSIDE OF PADDED SEAT!

ALSO AVAILABLE WITHOUT RETAIL PACKAGING. SEE PAGE 41

WEIGHS ONLY 12 LBS.

OPTIONAL BASKET AVAILABLE MODEL # 9501S72606

Aluminum Rollator with Fold Up and Removable Back Support, Padded Seat, 6" Casters with Loop Locks

- R726BL Blue, (Has replaced R6BL and 726BL), 1/cs
- R726RD Red, (Has replaced R6RD and 726RD), 1/cs
- R726SL Silver, (Has replaced R6SL and 726SL), 1/cs
- R726GR Green, 1/cs
- R726BK Black, 1/cs

FEATURES:

- 6" black, non marring casters with soft-grip tires are ideal for indoor and outdoor use.
- Comes with new seamless padded seat (14" x 14") with zippered pouch under seat.
- Removable, hinged, padded backrest can be folded up or down.
- Comes standard with carry pouch.
- Easy to use loop locks.
- Ergonomic handles are easy to grip and are adjustable in height.
- Comes with padded seat with zippered pouch under seat.
- Handles are adjustable in height.
- Limited Lifetime Warranty.

FOLD-UP & REMOVABLE BACKREST!



EASY HEIGHT ADJUSTMENT WITH SELF THREADING KNOB!

CERATED BRAKES



R726BL

PADDED SEAT OPENS TO A ROOMY CONVENIENT STORAGE POUCH! ALSO COMES WITH INSIDE ZIPPER POUCH ON UNDERSIDE OF PADDED SEAT!



COMES WITH NEW SEAMLESS PADDED SEAT WITH CARRY POUCH!



Blue

Red

Silver

Green

SPECIFICATIONS

SEE PAGE 52 & 53 FOR REPLACEMENT SEATS.

Item #	Handle Height (highest)	Handle Height (lowest)	Unit Width	Unit Length	Seat Depth	Seat Width	Seat Height	Product Weight	Carton Shipping Weight	Weight Limit
728RD-RTL	38"	33"	24"	28"	14"	14"	23"	14 lbs.	19 lbs.	300 lbs.
R726 SERIES	37"	32"	24"	28"	14"	14"	21"	12 lbs.	16 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.