



LEVER LOOP LOCK PROTECTS FINGERS FROM SNAPBACK!

EASY HEIGHT ADJUSTMENT WITH SELF THREADING KNOB!

WEIGHS ONLY 12 LBS.



PADDED SEAT MODELS COME WITH STANDARD POUCH (ITEM # 804)


SEE PAGE 52 & 53 FOR REPLACEMENT SEATS.

Deluxe Aluminum Rollator with Padded and Plastic Seat, 6" Casters with Lever Locks

- LBBL Blue, Plastic Seat with Basket, 1/cs
- LBRD Red, Plastic Seat with Basket, 1/cs
- LBPSBL Blue, Padded Seat and Pouch, 1/cs
- LBPSRD Red, Padded Seat and Pouch, 1/cs

FEATURES:

- Special lever lock prevents a person's fingers from "snapback" contact when releasing the wheel locks.
- 6" casters with soft-grip tires are ideal for indoor and outdoor use.
- Tool Free removable padded backrest with cover for comfort.
- Plastic seat models come with Large 14" (W) x 13" (D) x 6" (H) front mounting basket.
- Padded seat models come standard with pouch (item #804). The pouch adds privacy and security to personal items. Plastic seat models come with basket (Item # 830).
- Manufactured with wrap-around welds providing additional strength and support.
- Handles are adjustable in height.
- Easy, one hand folding on plastic seat models.
- Limited Lifetime Warranty.



PARTS DIAGRAMS ON PAGE 60

SPECIFICATIONS

Item #	Handle Height (highest)	Handle Height (lowest)	Unit Width	Unit Length	Seat (depth)	Seat (width)	Seat (height)	Product Weight	Carton Shipping Weight	Weight Limit
LB SERIES (PADDED SEAT)	38"	32"	23"	27"	13"	14"	21"	12 lbs.	16 lbs.	300 lbs.
LB SERIES (PLASTIC SEAT)	38"	32"	23"	27"	13"	14"	21"	12 lbs.	16 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.