

VIPER PLUS GT



EASY,
TOOL FREE
HEIGHT
ADJUSTMENT
ARMS



PADDLE ARM
RELEASE



QUICK RELEASE
AXLE

COMPARE TO SUNRISE® AND MEDLINE®

Aluminum Viper Plus GT Deluxe High Strength, Lightweight, Dual Axle (K0004)

FEATURES:

- Aluminum frame.
- Adjustable height, flip back, padded arms standard.
- VIPER PLUS GT weighs under 32 lbs. (Excluding front rigging). Easier to propel and transport than a standard weight chair.
- Adjustable angle back easily adjusts from 5° - 20°.
- Back height adjusts from 17" - 19" in 1" increments, providing added patient comfort and versatility.
- New frame style eliminates seat guides and allows for custom backs and accessories.
- Durable, heavy gauge flame retardant nylon upholstery with velcro modesty flaps, resists mildew and bacteria.
- **Quick release, 24" composite Mag-style wheels are lightweight and maintenance free.**
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability.
- Front caster forks have 3 height adjustments and angle adjustability.
- **Standard 3 position quick release axle provides easy transition of seat height.**
- Composite footplates are durable, lightweight and tool free height adjustable.
- Urethane rear tires and casters offer superior performance, smoother ride and are lightweight.
- **Adjustable anti-tipper with wheels flip up. Never needs to be removed.**
- Comes with push to lock wheel locks.
- Limited Lifetime Warranty.

NOW
COMES WITH
RIGGINGS IN
CARTON!

Due to product improvements and changes, all images in this catalog are subject to change.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.